# *Trim@gx*

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## Contents

---

**Caution:** Please read the following pages attentively.

---

|             |                                           |    |
|-------------|-------------------------------------------|----|
| Information | Indications .....                         | 25 |
|             | Contraindications . ....                  | 25 |
|             | Particulars of the device .....           | 25 |
|             | Warning and safety instructions .....     | 26 |
|             | Risks and side effects .....              | 27 |
|             | Maintenance and repair instructions ..... | 27 |
|             | Cleaning and care instructions .....      | 27 |
|             | Transport, storing and disposal .....     | 27 |

|                  |                                                  |              |
|------------------|--------------------------------------------------|--------------|
| <b>Treatment</b> | <b>Preparation of the device .....</b>           | <b>28</b>    |
|                  | <b>Preparation of the patient .....</b>          | <b>29</b>    |
|                  | <b>Operating the device .....</b>                | <b>29-33</b> |
|                  | <b>Positioning of the annular coils .....</b>    | <b>34-35</b> |
|                  | <b>Modulation of frequency .....</b>             | <b>36</b>    |
|                  | <b>Magnetising force at a coil diagram .....</b> | <b>37</b>    |

---

|                         |                                                  |              |
|-------------------------|--------------------------------------------------|--------------|
| <b>more information</b> | <b>Signs and symbols .....</b>                   | <b>38</b>    |
|                         | <b>Behaviour on malfunction .....</b>            | <b>39</b>    |
|                         | <b>Technical data .....</b>                      | <b>40</b>    |
|                         | <b>Warranty and liability restrictions .....</b> | <b>41-42</b> |
|                         | <b>Own notes .....</b>                           | <b>43-44</b> |

## Indications

---

The device can be used to support the pain reduction at arthropathies as well as to support wound healing.

## Contraindications

---

- Implanted active medical products as for example
  - Cardiac pacemaker
  - Heart defibrillators
  - Insulin pumps
  - Hearing aids
- Metal implants within the area of treatment
- Systemic infection
- Untreated Diabetes mellitus
- Bleeding duodenal and stomach ulcers
- Malignant tumours
- Pregnancy
- Epilepsy

## Particulars of the device

---

The device is controlled via touch screen. The integrated user's guidance leads through the menu step by step.


In the main menu two indication specific and one individual menu items are available to the user:

- Pain reduction on arthropathies
- Wound healing
- Manual operation

In the indication specific scope of applications 10 therapeutical sessions are stored, which differ in each case in the bandwidth of the frequency modulation.

Furthermore in the menu item >Manual Operation< the user has got the alternative to set the frequency and respectively the bandwidth of the frequency individually himself.

The duration of the therapeutical session is preset by the device to 20 minutes. This can be changed individually, if required.

The electromagnetic energy is applied to the area of treatment via annular coils. The annular coil is positioned in that way, that the aching joint respectively the wound is in the center of the ring. To reach each area of treatment optimally, there are 3 different annular coils available. Please see also "positioning of the annular coils".

The three annular coils are placed inside the case of the device. After choosing the needed ring size the user pulls the coil out of the case as far as required. After ending the treatment the cable can be coiled up again by pressing the key on the touchscreen and the coil can be placed again into the case of the device.

During this procedure the device can be locked by activating the brakes of the front rolls.

## Warning- and safety instructions

---

- The Trimagx may only be used for the kind of treatment mentioned in the operating instruction.
- The device may not be used in rooms, where flammable mixtures as there are anaesthetic agents, oxygen or nitrous oxid are located.
- The device may only be operated by persons qualified by their educational background and also instructed in detail in the handling of the device.
- The magnetic field therapy may only be used under survey of briefed personnel qualified in the adjuvant therapy of medical diagnosed diseases and disorders considering the contraindications
- Furthermore the contraindications must also be observed by the staff as well as third persons. A minimum distance to the active annular coils of 1 meter must be kept.
- The magnetic field therapy must not be administered simultaneously with other medical products (for example microwave therapy, electrotherapy, infusions).
- Please observe the instructions on possible risks and side effects.
- Please ensure that the annular coil is placed correctly at the area of treatment during the therapy.
- Electronic instruments and other medical products in the practice rooms, for example ECG-devices, can be interfered in their measuring and control function, if they are used at the same time. A minimum distance of 1 meter must be observed and the simultaneous operation of the devices must be avoided.
- The annular coils must be cleaned after each use.
- Receiver, for example radio-controlled watches or mobile phones, may be damaged by the influence of the electromagnetic energy. Please observe a minimum distance of 1 meter (see also chapter signs and symbols).
- During the treatment no pieces of jewellery may be worn by the patient in the area of treatment.
- Data media, as for example magnetic stripes on credit cards, can be deleted by the influence of the electromagnetic energy. A minimum distance of 1 meter must also be observed.
- Please see chapter signs and symbols
- In emergencies or at unusual reactions of the device the power supply plug should be pulled. After taking care for the patient please inform the after-sales service immediately.
- As supplier of the device we only can be held responsible for effects on safety, reliability and efficiency of the device, if
  - the device is only used solely according to the directives specified in our operating instructions
  - the electric installation of the relevant room meets with the DIN IEC/VDE standards.
  - The initial startup, maintenance, re-adjustments, modifications and extensions or repairs are only done by persons authorized by our company.
  - The medicine product book (according to MPBetreibV §7) is to be filled in and kept in due order.

## Risks and side effect

---

To avoid possible risks please observe the warning- and safety instructions as well as the before mentioned contraindications.

The following typical side effects can be noticed in individual cases:

- prickling within the area of treatment
- sensation of heat in the area of treatment
- alterations of sensibility
- sleep disorders
- light flashes

## Maintenance and repair instructions

---

**Technical safety controls may only be made by the manufacturer or other persons specially authorized by him.**

Because of the decree for the use of medical products (MPBetreibV) valid in Germany the user has to make sure that all necessary controls for functional and operating safety are effected (**§ 6.1**). The controls have to be made annually on the magnetic field therapy device. This control includes a visual inspection of the device and of the device specific components as well as a measurement of the electric safety acc. to VDE 701. This control may only be made by the manufacturer or persons authorized by him.

Also repairs may only be done by the manufacturer or persons authorized by him. On possible repairs please inform the manufacturer on mentioning the serial number of the device (type label on the back side of the device).

## Cleaning and care instructions

---

After each treatment please disinfect the annular coils with a moist cloth and the disinfectant Melisepto.

**CAUTION:** No liquids may enter the ventilation holes of the coils.

Furthermore the device and the annular coils may be cleaned only with water and a soft, moist cloth. Do not use any scouring agents or caustic cleaning agents.

## Transport, storing and disposal

---

The magnetic field therapy device may only be transported and stored in its original packing and at the temperatures indicated in the technical data sheet.

Old appliances can be disposed within the valid regulations on old appliances.


## Preparation of the device

Please observe the warning and safety instructions implicitly on choosing the location of the device.

The handle on the back side of the Trimagx has to be used together with the rolls to position the device. The handle is not suitable for carrying the device.

### !!! Caution !!!

Please remove the transportation lock before the first set-up (3x red screws)

The attached power cable is to be plugged in at the back of the device and connected to the power supply of the room. The device works with a line voltage of 230 V/ 50 Hz and may only be operated with the attached power cable. It is essential to activate the breaks at the rolls before pulling out the annular coils.


Turn on the device by using the switch at the back of the Trimagx.

The device will now carry out a self test.

ENGLISH

Afterwards the main menu appears on the touchscreen

## Preparation of the patient

---

Please bring the patient in a comfortable sitting or lying position, which enables an easy access to the area of treatment.

After choosing the needed ring size, please pull the coil out of the case as far as required.  
Above the touchscreen the symbol for the choosen ring lights up now.

S = small = 20 cm Ø

M = medium = 35 cm Ø

L = large = 55 cm Ø

The annular coil is positioned in that way, that the aching joint respectively the wound is in the center of the ring. Some examples are shown in chapter "positioning of the annular coils"


**Please take care that the ventilation holes of the coils won't be covered.**

## Operating the device


---

**On choosing an indication specific preassigned course of treatment:**

- Query on warning and safety instructions


- If you confirm this query with „no“ you will enter the warning and safety instructions.


- On each single button you get access to the respective information, on pressing the button „back“ you will access the main menu again.

- If you confirm the query on the warning and safety instructions with „yes“ you will enter the choice of the therapeutic session. You can change the quantity on activating the arrow keys

- If you have confirmed with „next“, you can now choose the duration of treatment. On using the arrow keys you can adjust the duration in steps of 5 minutes.

- If you have confirmed with „next“, the therapy data will be shown to you once again. Now you can pull out the suitable annular coil and start the treatment on pressing the button „Start the therapy“.


- During the whole therapy the current parameter will be indicated. You can stop the therapy at any time by using the button „Finish the therapy”.


#### On choosing manually an individual course of treatment :

- Query on warning and safety instructions


- If you confirm this query with „no“ you will enter the warning and safety instructions.


- On each single button you get access to the respective information, on pressing the button „back“ you will access the main menu again.

- If you have confirmed the query on the warning and safety instructions with „yes“, you can start the input of the individual frequency bandwidth. With the button “next” you will enter the frequency input.

- Determine the starting and final frequency on pressing the respective digit button.

- Via the number entry you can feed in the frequency and confirm it afterwards on using the enter key. With the key CLR you have the possibility to correct the entry.

- If you have fed in the frequencies and confirmed them with the enter key the therapy data will be shown to you once again. Now you can pull out the suitable annular coil and start the treatment on pressing the button „Start the therapy“.

- During the whole therapy the current parameter will be indicated. You can stop the therapy at any time by using the button „Finish the therapy“.


- Coiling in the cable and the ring will be done by activating the button „rolling in“.

## Positioning of the annular coils



---

- The annular coils are positioned in that way, that the aching joint respectively the wound is in the center of the ring.
- Therefore 3 annular coils in different sizes (Ø 55 / 35 / 20 cm) are available.
- During the treatment the patient should sit or lie as comfortable as possible.



### Examples:


Knee


Spine


Lumbar vertebra  
and hip




**Cervical spine**


**Shoulder**

**Shoulder and  
thoracic  
vertebra**

**Wounds**  
  
**at the foot**

**at the shinbone**

**ENGLISH**

**CAUTION:** Open wounds have to be covered sterily to avoid a direct contact with the coil.

## Modulation of frequency

The pulse frequency is one of the essential physical properties for the therapeutical efficiency of the magnetic field.

In general you can attribute a rather „calming“ or „sedative“ effect to the lower frequencies, while higher frequencies are more „stimulating“.

Furthermore you have to mind an intra- and interindividual variability concerning the effective frequencies during the treatment. Therefore the frequencies will be modulated in a specific bandwidth at the identification specific preassigned programs.

**Example:** The frequency value 100 Hz means, that the magnetic field is build up and relieved 100 times per second.

On choosing a frequency bandwidth of 100 - 200 Hz the therapy program starts with the value „100 times per second“ and will get „faster“ during the therapy session: 101 Hz, 102 Hz, 103 Hz and so on.

As soon as the program has reached 200 Hz the frequency will slow up until it has reached the initial value of 100 Hz.

This procedure will be repeated until the end of the therapy session.


### Frequency control scheme:

| Arthropathy   |                  |
|---------------|------------------|
| Treatment No. | Modulation<br>Hz |
| 1             | 1 – 100 – 1      |
| 2             | 1 – 200 – 1      |
| 3             | 1 – 50 – 1       |
| 4             | 1 – 75 – 1       |
| 5             | 1 – 200 – 1      |
| 6             | 1 – 200 – 1      |
| 7             | 1 – 75 – 1       |
| 8             | 1 – 50 – 1       |
| 9             | 1 – 200 – 1      |
| 10            | 1 – 100 – 1      |

| Wounds        |                  |
|---------------|------------------|
| Treatment No. | Modulation<br>Hz |
| 1             | 1 – 5 – 1        |
| 2             | 1 – 15 – 1       |
| 3             | 1 – 75 – 1       |
| 4             | 1 – 5 – 1        |
| 5             | 1 – 75 – 1       |
| 6             | 1 – 100 – 1      |
| 7             | 1 – 5 – 1        |
| 8             | 1 – 100 – 1      |
| 9             | 1 – 75 – 1       |
| 10            | 1 – 100 – 1      |

## Magnetising force at a coil diagram

Coil diagram – large coil diameter 55 cm

Each colour ring marks a magnetising force.

Initial point is the grey ring, which has a force of 30 mT at each of the 3 annular coils large/medium/small.

The magnetising force decreases constantly inwards and outwards.

- |     |               |                          |
|-----|---------------|--------------------------|
| 7.  | Zone = grey   | 30 mT (at the coil body) |
| 8.  | Zone = orange | 25 mT                    |
| 9.  | Zone = blue   | 18 mT                    |
| 10. | Zone = green  | 10 mT                    |
| 11. | Zone = yellow | 4 Mt                     |
| 12. | Zone = lilac  | 0,4 mT                   |

## Signs and symbols

At the upper side of the device there are 3 indicator lamps above the touchscreen. They have the following meaning:

**S** ring coil „small“ removed  
**M** ring coil „medium“ removed  
**L** ring coil „large“ removed

Observe operating instructions

**Sicherheitsabstand 1 Meter**  
safety clearance 1 meter


Safety and after-sales service

Warning of hand injuries

No metallic parts or watches allowed


Please observe the legal prescriptions (see maintenance instructions)

Prohibition for persons with cardiac pacemakers


No mobile radio allowed

Disposal (pls. see transport/storing/disposal)

Prohibition for persons with metallic implants

CE-sign of the mentioned board

No magnetic or electronic data media allowed

## Behaviour on malfunction

---

Should any malfunction occur while starting or operating the device you will get the following information at the touchscreen:

Please follow the described procedure:

The touchscreen is equipped with a standby-function. If the device won't be used in switched-on condition for 30 minutes, the touchscreen will switch off automatically. For reactivating it please touch the surface of the touchscreen slightly.

## Technical data

---

|                                                      |                                                              |
|------------------------------------------------------|--------------------------------------------------------------|
| Name of appliance:                                   | Trimagx                                                      |
| Type of appliance:                                   | Magnetic field therapy device<br>UMDNS-Nr. 12415             |
| Power supply:                                        | 230 V~/50Hz                                                  |
| Power input:                                         | 470mA                                                        |
| Fuses:                                               | 2 x T1,6 A ,L 250V                                           |
| Classification:                                      | Schutzklasse I                                               |
| Dimensions of device:                                | 535 x 1130 x 770 mm                                          |
| Dimensions annular coils:                            | 550 mm, 350 mm, 200 mm                                       |
| Maximum power output:                                | 30 mT at the coil body                                       |
| Connections:                                         | 3 annular coils with size S/M/L                              |
| Power line:                                          | Item No. 40126                                               |
| Ambient temperature                                  |                                                              |
| a) during operation:                                 | +15 °C - +35 °C                                              |
| b) at stock:                                         | + 5 °C - +50 °C, relative humidity < 90% (not<br>condensing) |
| Atmospheric pressure:                                | 600 – 1200hPa                                                |
| Weight:                                              | 78 kg                                                        |
| Product category acc. to<br>MP-directives 93/42 EWG: | II a                                                         |

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## Warranty conditions

---

The invoice, as a receipt for your purchase, should be kept carefully. You will need it for warranty claims.

HERRMANN Apparatebau GmbH guarantees that the product is fully functional and technically concordant with the documentation.

During production and before delivery, our products undergo high quality checks. If there is yet a deficiency or processing fault, HERRMANN Apparatebau GmbH issues a warranty of two years, valid from date of first purchase. During this term HERRMANN Apparatebau GmbH remedies deficiencies gratis, which demonstrably result from material- or processing faults either with repair of the product or replacement of faulty parts.

### The warranty does not apply:

- with deficiencies which result from transportation, accidents, natural disaster, vandalism, improper use, incorrect maintenance or incorrect repair through others.
- with change of software (ozon devices), maloperation, incorrect installations, or other modifications not agreed with us.
- with disrespect of instructions in enclosed documentations.
- with incompatibility or malfunction, resulting from not through us used product components.
- with phenomena which occur with the regular aging process of the product.
- with deficiencies resulting from external appliances.
- with any deficiencies of a product whose serial number has been changed or removed.

Warranty terms for replaced parts end together with the end of the original warranty term for the complete product.

If you return a device on terms of warranty, products sent to us without any accessories will be replaced without any accessories. If you return the product without its original packing and the respective packing cushions HERRMANN Apparatebau GmbH is not liable for all damages, which are accrued during transport. HAB will only be liable for the complained deficiencies or faults

If you assign your device to a different user, the remaining term will be assigned to the new user. Invoice, warranty operating instruction and medical device book should accompany the device.

## Liability restrictions

---

Claims for damages are excluded unless they are based on intention or gross negligence of HERRMANN Apparatebau GmbH or its employees. The liability according to the law for product liability will remain unaffected. On no account HERRMANN Apparatebau GmbH is liable for:

- Claims raised against you by third persons due to loss or damage.
- Economical consequential damages (including lost winnings or economies) or accompanying damages. That also applies to cases in which HERRMANN Apparatebau GmbH was informed about the possibilities of such damages.

HERRMANN Apparatebau GmbH is on no account responsible for any random, indirect, special consequential or other damages (including, without any restrictions, damages in terms of loss of profit, business interruption, loss of business information or any other losses), which arise with the use of the product or in any relation to the product based on a contract, compensation, negligence, strict liability or other claims also if HERRMANN Apparatebau GmbH was informed in advance of the possibility of such damages. This exclusion also includes every liability, which can result of third persons claims against the first buyer. This sale clause aims to restrict possible liability of HERRMANN Apparatebau GmbH against claims which can result from the restricted warranty and/or the sale.

The above declaration does not take effect in some countries where this exclusion or the restriction of accompanying or consequential damages is legally not allowed.

## Own Notices

---